




mCon1/JFHrluc2(j5.1; adapted)

Virus Specifications:

Description:

mCon1/JFHrluc2(j5.1; adapted) is a genotype 1b/2a chimeric reporter virus in which the HCV IRES drives translation of a *Renilla* luciferase-containing HCV polyprotein that is cleaved by cellular and viral proteases to produce the individual HCV proteins. The *Renilla* luciferase is fused C-terminal to the p7-NS2 cleavage site, and is cleaved from the HCV polyprotein by virtue its fusion with the foot-and-mouth disease virus 2a peptide coding sequence (FMDV2a). This virus contains five adaptive mutations that were identified after passaging the parental virus in cell culture to generate a variant that produces high titers.



Genotype: 1b/2a

Strains: Con1—core through the N-terminal region of NS2
JFH1--5'NTR; N-terminal region of NS2 to NS5B
J6--3'NTR¹

Location of chimeric junction: between the first and second predicted transmembrane domains of NS2²

Reporter: *Renilla* luciferase

Configuration: monocistronic

Reporter location: C-terminal to the p7-NS2 junction

Adaptive mutations: one within the Con1 region of NS2, one within the JFH1 region of NS2, one within NS3, two within NS5A³; adaptive mutations are indicated in the above figure with an asterisk.

Maximal titer range: low 10³ to low 10⁴ TCID₅₀/ml

Time point post-electroporation at which maximal titers are achieved: 48 hours

Plasmid designation: APP152

Virus designation: APV152

Virus Production:

Description: Due to the error-prone nature of the HCV RNA-dependent RNA polymerase, it is recommended that stocks be generated from *in vitro* transcribed (IVT) RNA generated from the cloned viral genome. Virus is produced following electroporation of IVT RNA into the highly permissive human hepatoma cell line Huh7.5. The Huh7.5 cell line was generated by curing a stably selected replicon-containing cell line⁴, and are the industry-standard for the production of high titer culture-produced HCV.

Huh7.5 cell line designation: APC49

Cell yield from a confluent T150 flask: 6 – 8 x 10⁶ cells

Growth conditions: 5% CO₂ and 37°C

% viability post-freeze: > 90%

¹ Differences between the J6 and JFH1 3'NTRs lie within the variable region; side-by-side comparison of a construct containing the J6 3'NTR versus the JFH1 3'NTR revealed no titer differences

² Pietschmann *et al*, *PNAS* (2006) 103: 7408-13

³ Lacks adaptive mutations in E1 and E2

⁴ Blight KJ *et al*, *J Virol* (2002), 76: 13001-14




Figure 1. Huh7.5 cell titration growth curve. The indicated number of Huh7.5 cells was seeded into 96-well plates at time = 0. To quantify the relative number of viable cells, intracellular ATP levels were quantified 4, 24, and 48 hours post-plating.

Recommended HCV Antiviral Assay:

1) Infection assay to determine efficacy (EC₅₀) and viability (CC₅₀):

EC₅₀ assay: Test Article (at varying concentrations) and HCV virus that expresses *Renilla* luciferase are added simultaneously to Huh7.5 human hepatoma cells. *Renilla* luciferase levels are quantified 48 hours after Test Article/virus addition using the Promega *Renilla* Luciferase Assay System to determine the level of viral replication inhibition.

CC₅₀ assay: Test Article (at varying concentrations) is added to Huh7.5 human hepatoma cells in the absence of HCV virus. Intracellular ATP levels are quantified 48 hours after Test Article addition using the Promega Cell Titer-Glo^R Luminescent Cell Viability Assay System.

Raw data format EC₅₀: *Renilla* luciferase expression in Light Counts Per Second (LCPS)

Raw data format CC₅₀: Luminescence measured in LCPS

Calculated data format: Percent Inhibition relative to solvent vehicle-treated control

Maximum number of EC₅₀ and CC₅₀ values per 96-well plate: four

Signal-to-noise ratio: >500

Efficacy Determination for HCV Inhibitors:




Figure 2. EC₅₀ determination in APV152 infection assay in which the Test Article is a protein-based HCV entry inhibitor. Cells were seeded 24 hours prior to addition of Test Article and APV152 at MOI = 0.1. The assay proceeded for 48 hours, at which time *Renilla* luciferase activity was measured. EC₅₀ = 0.12 ug/ml




Figure 3. EC₅₀ determination in APV152 infection assay in which the Test Article is an HCV polymerase inhibitor. The assay was performed as described in Figure 2. EC₅₀ = 0.12 uM, CC₅₀ >25 uM (data not shown)




Figure 4. EC₅₀ determination in APV152 infection assay in which the Test Article is an HCV protease inhibitor. The assay was performed as described in Figure 2. EC₅₀ = 0.28 uM, CC₅₀ > 25 uM (data not shown)




Figure 5. EC₅₀ determination in APV152 infection assay in which the Test Article is an HCV cell-based inhibitor. The assay was performed as described in Figure 2. EC₅₀ = 0.15 uM, CC₅₀ > 25 uM (data not shown)